

CCDC43-coiled-coil domain containing 43 Gene

Reference: PTXBC047776

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coiled-coil domain-containing protein 43; coiled-coil domain containing 43		
For research use only.			

GENE DESCRIPTION	coiled-coil domain containing 43
GENE ID	124808
ORIGIN SPECIES	Human
ACCESSIONS	BC047776
VECTOR	pDONR223
NCBI SYMBOL	CCDC43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgccagcgaagtggcgcgatagccctggcgaagcgcgatggcggcggcgtttggctcctggcggcggcgttggaggcactggagtgaccgagccgtttatggagcctacatctgggtatcctgcaggaggaggaggagaagagaagctggac
gctctgcaggggatcctctctgtttcttggagaagatccctccttaatatctgcaaggagattgtggaacgatggtcagaaactcagaatgtgtcacaaagtgaagaaagagatgaagtacaggccattgccacctaattgagaagcaggcacaatctagttaagccaagatg
gtgtcagaagaggagaagcagagaaaaagctgccctcctggccagatgctgacagatgaaggaggaatcacttactgaagcagatgagaaggatgattcaggtgtctaccacaatgaacattggttctgacaaactctgttccgaaacccaatgtggaagatgcctaatgccg
aaaactggagcgagactcctcgggatgaatccaaaggagaaggaacaggacaagctgcagagggaagagacaaactagccaagcaggagcgaaggaaaaggaaaaaaaggacacagagagggagcgaagcgataa