

ARHGAP24-Rho GTPase activating protein 24 Gene

Reference: PTXBC047918

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FILGAP; RC-GAP72; RCGAP72; p73; p73RhoGAP; rho GTPase-activating protein 24; RAC1- and CDC42-specific GTPase-activating protein of 72 kDa; filamin-A-associated RhoGAP; rhoGAP of 73 kDa; sarcoma antigen NY-SAR-88; Rho GTPase activating protein 24		

For research use only.

GENE DESCRIPTION	Rho GTPase activating protein 24
GENE ID	83478
ORIGIN SPECIES	Human
ACCESSIONS	BC047918
VECTOR	pDONR223
NCBI SYMBOL	ARHGAP24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagaacaatgactccacggagaaccccccaagccaaggcgccagaaatgccatcaagtgtgggtggctgagggaagcaaggagccttgtcaagactggcactacgcgtgtgtgtctcaagggggatcagctctattattcaaatgaagatgaacaaagcccttg
ggfataatttttgcctggaataaagtgtctgagcatccctgcaatgaagagaacccagggaagttccttttgaagtagttccaggagcgcgatcgagatcgatgacagcaaatcatgagagctacctctcatggcaagcaccagaatgatafgaagactgggtgaagtcaatccgc
cgagtcataatgggacctttcgaggaggcatttttggacagaaactggagatactgttcgtatgagaagagatatgggaaccgtctgctccgatgtgtgtggagcagtcgctggacttatccgacaaaggggctgaagaagagggtctcttctgactgccaggccaggcctaactc
tgttaaggagctccaagatgcccttgactgtgggagaagccatcattgacagcaacacagatgtacacacgtggcatcacttctaagctgtacctccgagaacttcagaaccagttattcctatcggaagtatgaagattttgtcatgtgccaaactgctcagcaaggaggaag
cagtaagttga