

## VPS4A-vacuolar protein sorting 4 homolog A (*S. cerevisiae*) Gene

Reference: PTXBC047932

|                |                                                                                                                                                                                                            |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                        | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                       | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | SKD1; SKD1A; SKD2; VPS4; VPS4-1; vacuolar protein sorting-associated protein 4A; SKD1-homolog; hVPS4; vacuolar protein sorting factor 4A; vacuolar sorting protein 4; vacuolar protein sorting 4 homolog A |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | vacuolar protein sorting 4 homolog A ( <i>S. cerevisiae</i> )                                    |
| <b>GENE ID</b>          | 27183                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC047932                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | VPS4A                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgacaacgtcaacctccagaaagccattgatctgtgtgacgaaagccagaggaggacaaagccaagaactacgaggaggcgctgcggctgtaccagcatcggtggagtactctccacgctatcaagatgaggccacagcgacaaggccaaggagagcattcgagccaa  
gtgcgtgcagtagaccgggcccagagaagctgaagatatttacgaagcaagagaaacacggcaagaagccagtcaagagaaaccagagtgaggggcaaggcgagtgacagtgacagtgaaagggaataatccggagaaaaaagaaactgcaagaacagctgatgggtgcctg  
cgtgatggagaagcccaacatagctgtggaacgacgtggccgggctggagggggccaaggaggccctcaagaagctgtcatctttgccaatcaaatccacactgtgtcacaggcaagcgacccccctggcggggattctgtctgttcggacccccctggcacagggaatcctacctg  
gccaaagccgtggcaacagaggccaactccactctctctgtgtcctcctcagatctgatgtccaaagtggctggggagagtgaaaagctgtgtcaagaacctgtttgagctggccaggcgacacaaagccctccatcttcacgatgagggtggaattccctctgcgggtccgaaa  
tgaaaatgagagtgaggccgcccggagatcaaaacggagttctgtgtccagatgcagggggtggggaataacaatgatggactctgttctgtggagccacaacatccatgggtgttggaattcgccatcaggaggaggttgaaaaacgaatttatcccttgcgggaggaagc  
tgcccgccagatgttccggttgcatctcgggagcactccccacaacctcaggatgcaaacatccacgagctggccgggaagacggaaggctactcggcgcgacatcagcatcatcgtgcgggactctctcatgcagccccgtgaggaagggtgcagtcggccacactcaa  
aaaggtctgtggccctctcgcacaaacccagcatgatgattgagactcctgactccatgctcaccaggggacccagagccatgagatgacttgatggatgtccctggggacaaactcttagagcctgtgtgttcagtgctggacatgctgcggtctctgcccaccacccggccc  
acgggtgaatgcagacgacctcctgaaagtgaagaatttcagaggactttggccaagagagttaa