

GRID1-glutamate receptor, ionotropic, delta 1 Gene

Reference: PTXBC047948

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GluD1; glutamate receptor ionotropic, delta-1; gluR delta-1 subunit; glutamate receptor delta-1 subunit; glutamate ionotropic receptor delta type subunit 1		

For research use only.

GENE DESCRIPTION	glutamate receptor, ionotropic, delta 1
GENE ID	2894
ORIGIN SPECIES	Human
ACCESSIONS	BC047948
VECTOR	pDONR223
NCBI SYMBOL	GRID1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacagcctcatggaagacattgctcacaagcagattccccagcgtcgattgagctctcgccctggagatggggggcctggctccaccagaccttgagccgacacgggagtagcagaacaccagctctcggtcagcacctttctccagagcagagcagccatggca
ccagccggacactctcatcagggccagcagcaacctgcccgtgcccgtgagcagctcgccgacctgccctccatgcagtgcacaacacaggtcacccaacggggggctgttcggcagagcccgggtgaagacccccatccccatgtccttcagcccgtgcttgaggcgcttc
cagaggctctggacacctccacgggacctccatctga