



Reference: PTXBC047949

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM102A; AI426465; C230093N12Rik; Eeig1; early estrogen-induced gene 1 protein; family with sequence similarity 102, member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 102, member A
GENE ID	399665
ORIGIN SPECIES	Human
ACCESSIONS	BC047949
VECTOR	pDONR223
NCBI SYMBOL	FAM102A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgtctctgctctctggaatccctgcttcaacagcgcacacatgactgcgaagtcacatctccatccaggccaggaatctctctctgctgactgacgtgtaaggtgtgtgggaccagcagtgggggcagcagcaccaactcctgactgggtccggcccccgaagctcggcccatatt
ctcagctcagggctgccaagggaacccgaccagaacctgtccagccctgaggagggtgttccactctggccactccgcaactcagctatgccagccagcagcttcaagatctccggctacagcacagagcactgcgctctccagctctcagacctgacgcaccgccgcaacagct
caccagcagcagcgctctgggggcttggcgtgacgtctgagggcccttgagggcagtgagcggagaccggcccccggagaagccggccggcgccaccgccggccctgcatctgtccatgctcttttcaggcggaagaaggaactgggtfagagaccacccgactgggtg
gtagacacgcggatgatgctggatgcatctgtggagaagatctgcagagccaggaactcacagatggcagcaaacccaggagacacactccggctgtctgtgagcccgatggctctgccacgtgagggcatccactgtgcccacgtgctcctctctgggtctacagcca
cttctgattgaagccattga