

ZSCAN1-zinc finger and SCAN domain containing 1 Gene

Reference: PTXBC048099

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MZF-1; ZNF915; zinc finger and SCAN domain-containing protein 1; zinc finger with SCAN domain 1; zinc finger and SCAN domain containing 1		

For research use only.

GENE DESCRIPTION	zinc finger and SCAN domain containing 1
GENE ID	284312
ORIGIN SPECIES	Human
ACCESSIONS	BC048099
VECTOR	pDONR223
NCBI SYMBOL	ZSCAN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttccacgcccgaagccctgctccccagacgccccagacccaaccccagtgagcaggacgcagaccctggccagcaagcccgaggacacgaagccagcgtctgcgttcggcagttccagttaccagtgaggagcgggcccacctggcgtgggcca
gctctggacgctgtgcccagtggtgagccagaggcgcgtccaaaggagcagatgctggagctgctgctgagcagttctggcgctgcccagcaagatcgggacctgggtgcagtcacaggccccgaagctgcaggaggcccccagcctgtggaggacct
cacacagatgtgccagcaggaaggtgagaggcgcaggcttcgtgccccggcgccagggggcctgacgtctccgaggaccgaactgcctctccgcccacatcgctccaccctgcccggccctcccgaccagcaaacctctctgccccggcgtgttttaattatt
tttacatataa