

ZNF333-zinc finger protein 333 Gene

Reference: PTXBC048107

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 333		
For research use only.			

GENE DESCRIPTION	zinc finger protein 333
GENE ID	84449
ORIGIN SPECIES	Human
ACCESSIONS	BC048107
VECTOR	pDONR223
NCBI SYMBOL	ZNF333
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaatccgtcaccttgaggatgtggccgtggaggtcatccaggagtgaggcattgtggacacgcacggaggagcctgtgcaaatacagatgcttgaccagtgcaggaccctggcctccagggggaactccaccatgcaaacccagttgtgtctccagctggggcaagagcag
agccaaaggcaacagaacgagggtatctccgtgcccacaggtgtgcttgggaatcctaactaaacccgaagaggtgtgcttctatgcaggatcttttgaagaagcatcctccaggacatgcaaatggggccggggctgttctgaggatgcagctgggtcccccatagaagagagg
gagacacattgactcagaggaccggccagctctccaggagccgcttggctctgggatgcacgggactgaaggccgctatgcagattcagaggggtgtgataccagtgcctactctgggccaccgcaacccatgggtggccagggtattctgtgtgctgcacgtgacctgtgctg
gcttcaggaggacaaagtggaggaagctatggtctctgggtgccaacgcctgtgtcacaggaaaccagtcacctttgcagatgtggctgtgtgttcacccagaagaatgggtgtttctggaactctactcagaggagcctgtatagagatgtgatgtggagaactacaggaaactg
gcctctgtgtgatcaactgtgcaaaccaatgcgtgtctatttggagaagaggagagcagtgaccactgacagggcgctcctctcagacactgtgcagaacctcagtgcaacccaagaggaatcctagccaagatactttacagagatcctgtcattgatgtgaagggg
ccaccaactaa