

## C15orf33-chromosome 15 open reading frame 33 Gene

Reference: PTXBC048125

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C15orf33; protein FAM227B; family with sequence similarity 227 member B		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	chromosome 15 open reading frame 33
<b>GENE ID</b>	196951
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC048125
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	C15orf33
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

```
atggcaggccaacgaacatgtcaaaaggaggagcagcagactggccccgggaaatgcaagaacctcaaaagacattgaagaattctaaagttcaaaatgggattatggccaagggaatccatttagagatgatataaatggtcagcacttgaataaaataaaagaata
gttcattgtttcaattatacacacctatgggaaatgttcctcgaatattgaagcactttgatcatggaatcaaaataaaaggaatattcacttattttcaaaaccacacctctgaatatttaagtggaaaagcatgattctgagacaagttctatagaaaatggaacgatatgggaaattcta
aagaagtaccataagaaaaaaagataatgctttcagatgaaatggagacagagaagaatagagggttcagcgttcacaggattcaaaacaaatgaacttactcagctaccagacattggatgctgaacaaattatcttttattttaaaagccataatttgatgaaagattttaaaat
ctggaagactcattttctcagaagcctccattgctcttfgcatgactcctttgggtgggtttctccataaattaggatatactgattgtggcacaagccatataccaacgttccatgaagcatttccagaatcagttacctcttfaatgatgaattaaagaagatctagggaataacattt
tctttggttcagggtttaaaacctcaaaaaggcttttgatccactgaaactgaaagaactctcacaaccaccattcatgtagcaaaaagcacctgcaaaaatcagtaaaaggaaagaattgcagacagcaagaacatatacaactagctga
```