



Reference: PTXBC048285

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Lyl10; TSCOT; thymic stromal cotransporter homolog; thymic stromal co-transporter; solute carrier family 46 member 2		

For research use only.

GENE DESCRIPTION	solute carrier family 46, member 2
GENE ID	57864
ORIGIN SPECIES	Human
ACCESSIONS	BC048285
VECTOR	pDONR223
NCBI SYMBOL	SLC46A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]