

SYS1-SYS1 Golgi-localized integral membrane protein homolog (S. cerevisiae) Gene
Reference: PTXBC048286

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SYS1, golgi trafficking protein; Sys1 golgi trafficking protein; SYS1 Golgi-localized integral membrane protein homolog; protein SYS1 homolog; C20orf169; dJ453C12.4; dJ453C12.4.1		

For research use only.

GENE DESCRIPTION	SYS1 Golgi-localized integral membrane protein homolog (S. cerevisiae)
GENE ID	90196
ORIGIN SPECIES	Human
ACCESSIONS	BC048286
VECTOR	pDONR223
NCBI SYMBOL	SYS1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggtcagttccgacgtacgtgtgggacccgctgctgatcctgtcgcagatcgtcctcatgcagaccgtgtattacggctcgtggcctgtggctggcgtgtgggacgggctagtgcgaagcagccctcgtggaccagatgttcgacgccgagatcctgggctttccaccc
ctccaggccggctctccatgatgtccttcacccaagccctcaactgtgcccgtggcctgtgtgtactcaccggcgaggaaagcagtgctggtttcactgtcactgtccatttcttcacctctgggctgctgtgttctacagctccgtttccctcggcgctgacctgggtggtggtccaa
gccgtgtcattgcactcatggctgtcatcggggagttacgtgtcatgcggacggagctcaaggagataccctcaactcagccctaaatcaatgtctag