

SGOL2-shugoshin-like 2 (S. pombe) Gene

Reference: PTXBC048349

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SGOL2; TRIPIN; shugoshin 2; shugoshin-like 2		
For research use only.			

GENE DESCRIPTION	shugoshin-like 2 (S. pombe)
GENE ID	151246
ORIGIN SPECIES	Human
ACCESSIONS	BC048349
VECTOR	pDONR223
NCBI SYMBOL	SGOL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtccccagtgatggaactggctcactttttacctcaggaattaagagacatttgaagacaaaagaatttcaagactactaagtgaatgtttctcttcttcaaaaaataaaacaaaaataactaataattctctattttcaaaatattttaagcacaacaacagggcattagctcagg
ctcttagtagagaaaaagagaattctcgaagaattacaactgaaaagatgctatttcaaaaaagaagtagagaaactgaattttgagaacacatttctcctaaagctaaataacttgaataagaagcttatagacatagaagctctcatgaacaataacttgataactgcaattgaaatgagca
gtctttcaggttccatcagagttccttctactgtcagctagcaagaagaacgaattagtaaacagtgcaggttgatgcgtcttccatttgcagggttccattcaaatgatgatgaagatgaagataaagagaaaaatgcagtgtgacaacaataataatcaaaacattacctgatattc
cctcttcaggatcaacaacaacctttatcaactcaggataattcggagtggttatttctaaagaaaaatacaaatgtatatgtgttagatgattcagaacataattctctatagtgtgatgtacctccagaggtgagattgttttaaaattcattttgagtatttgattag