

C17orf64-chromosome 17 open reading frame 64 Gene

Reference: PTXBC048806

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C17orf64; chromosome 17 open reading frame 64		
For research use only.			

GENE DESCRIPTION	chromosome 17 open reading frame 64
GENE ID	124773
ORIGIN SPECIES	Human
ACCESSIONS	BC048806
VECTOR	pDONR223
NCBI SYMBOL	C17orf64
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctctggcgggtcatctccctctctcggagctggaagcaaaagcagcttcgccgactttacaagtacacaaagcagccagccgccaagttcctgtgacattctgcgcctcggatgcgcggagaggtccttgcctggccgaccgggaagacagcttgcctcaagctctgccatgcat
gggggctgcacagcaatcagcggcatgaaggagcggctgtccaatgcagacccaggtcaaggaggccccctgcctggcagccaagatcccaggacatgtcaagaaagattcttaaggagctttctcaaaaacaaaactcaagaggaaggataaaggagcccca
gaaactccagagactgaaccgtaa