

INO80C-INO80 complex subunit C Gene

Reference: PTXBC048989

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C18orf37; IES6; hIes6; INO80 complex subunit C; IES6 homolog		

For research use only.

GENE DESCRIPTION	INO80 complex subunit C
GENE ID	125476
ORIGIN SPECIES	Human
ACCESSIONS	BC048989
VECTOR	pDONR223
NCBI SYMBOL	INO80C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagccatgagtgaataaaatggcctctgagtttagcacaggacctgtggaaaaagctgccaaacctttgccatttaaggatcccaactttgtgcactctggccacgggtggcgcagtagctggcaagaagaacagaacctggagaacctgaaacaaatcctcgttctgaaa
gggcattgccgtggcaactgaacgatcctaactactcagtaattgatgctcctccatcctttaagccagctaagaagtattctgatgtttcagggtctgtgccaaactacacagacccccagagcaaacctgcggttcagcaccattgaagagttttctacattcggaggctgccctctgacgtcgt
caccggctacctggcctgagggaaggccacgagcatcgttcctga