



Reference: PTXBC048991

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	6M1-16; HS6M1-16; OLFR42A-9004-14; OLFR42A-9004.14/9026.2; OR2H6; OR2H8; OR6-2; dJ994E9.4; olfactory receptor 2H1; olfactory receptor 2H6; olfactory receptor 2H8; olfactory receptor 6-2; olfactory receptor OR6-32; olfactory receptor, family 2, subfamily H		

For research use only.

GENE DESCRIPTION	olfactory receptor, family 2, subfamily H, member 1
GENE ID	26716
ORIGIN SPECIES	Human
ACCESSIONS	BC048991
VECTOR	pDONR223
NCBI SYMBOL	OR2H1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]