



Antibodies and Proteins

## ProteoGenix SAS

## Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

[www.ProteoGenix.science](http://www.ProteoGenix.science)

## USF2-upstream transcription factor 2, c-fos interacting Gene

Reference: PTXBC049821

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

**Delivery Condition:**

+4°C, Blue Ice

### Synonyms

FIP; bHLHb12; upstream stimulatory factor 2; c-fos interacting protein; class B basic helix-loop-helix protein 12; major late transcription factor 2; upstream transcription factor 2, c-fos interacting

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | upstream transcription factor 2, c-fos interacting                                                |
| <b>GENE ID</b>          | 7392                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC049821                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | USF2                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

[illegible]