

CD3E-CD3e molecule, epsilon (CD3-TCR complex) Gene

Reference: PTXBC049847

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD3e molecule; CD3e molecule, epsilon (CD3-TCR complex); CD3e antigen, epsilon polypeptide (TiT3 complex); IMD18; T3E; TCRE; T-cell surface glycoprotein CD3 epsilon chain; CD3-epsilon; T-cell antigen receptor complex, epsilon subunit of T3; T-cell surface		

For research use only.

GENE DESCRIPTION	CD3e molecule, epsilon (CD3-TCR complex)
GENE ID	916
ORIGIN SPECIES	Human
ACCESSIONS	BC049847
VECTOR	pDONR223
NCBI SYMBOL	CD3E
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagtcgggcactcactggagagttctggcctctgcctcttaccagttggtgtttggggcgaagatggaatgaagaaatgggtggtattacacagacaccatataaagttccatctctggaaccacagtaaatatgacatgccctcagttatcctggaataactatggcaacacaat
gataaaacatagggcgtgatgaggaatgataaaacatagggcagtgatgagatcacctgtcactgaaggaaatttcagaattggagcaaaagtgggtattatgtctgtacccagagggaagcaaacagaaagatgcgaactttatctctactgagggaagagtggtgagaactgcatg
gagatggatgtgatgctgggtgccacaattgcatagtgacatctgcatcactggggcctgctgctgctgtttactactggagcaagaatgaaaggccaagcctgtgacacgaggagcgggtgctgctggcaggcgaaggggacaaaacaaggagagccaccact
gttccaaccagactatgagccatccggaaggccagcgggacctgtattctgcctgaatcagagacgcatctga