



Reference: PTXBC050312

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AKAP-13; AKAP-Lbc; ARHGEF13; BRX; HA-3; Ht31; LBC; PRKA13; PROTO-LB; PROTO-LBC; c-lbc; p47; A-kinase anchor protein 13; A kinase (PRKA) anchor protein 13; A-kinase anchoring protein; LBC oncogene; breast cancer nuclear receptor-binding auxiliary protein;		

For research use only.

GENE DESCRIPTION	A kinase (PRKA) anchor protein 13
GENE ID	11214
ORIGIN SPECIES	Human
ACCESSIONS	BC050312
VECTOR	pDONR223
NCBI SYMBOL	AKAP13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]