



Reference: PTXBC050387

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ankyrin repeat and SAM domain-containing protein 3; ankyrin repeat and sterile alpha motif domain containing 3		

For research use only.

GENE DESCRIPTION	ankyrin repeat and sterile alpha motif domain containing 3
GENE ID	124401
ORIGIN SPECIES	Human
ACCESSIONS	BC050387
VECTOR	pDONR223
NCBI SYMBOL	ANKS3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]