

HOXC4-homeobox C4 Gene

Reference: PTXBC050442

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HOX3; HOX3E; cp19; homeobox protein Hox-C4; homeo box 3E; homeobox protein CP19; homeobox protein Hox-3E; homeobox C4		

For research use only.

GENE DESCRIPTION	homeobox C4
GENE ID	3221
ORIGIN SPECIES	Human
ACCESSIONS	BC050442
VECTOR	pDONR223
NCBI SYMBOL	HOXC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcatgagctcgtatttgatggacttaactacatcgatccgaatttctccatgcgaagaatttcgcaaaatagctacatccctgaacacagtcgggaattacggccggaccagggaatcgggattccagcatcaccacaggagctgtaccaccaccgctccgcccctagc
taccctgagcggccagtatagctgcacagctccaggggccggcaattcgcgagggccacggccggccagggccggccaccacaccccgagaatcacagctcgtctcgagccggcgctctctcaggcgctcgcctccccgctccccagccccgccagctgcagccagc
cagccccgaccatccctccagcggccagcaagcaaccatagctacccatggatgaaaaaaatcacgttagcacgggtgaacccagttataacggagggaacccaagcgtcagggacagcctatacccgagcgaagctcgtgaattagagaagagtttcaatacaaccg
ctacctgaccggaaggagaaggtcagatcgccactcgtgtgctctctgagagggcagatcaaaatctgttccaaaacgctgcgatgaaatggaagaaggaccacgactcccaacaccaaaagtcaggtcagcaccgccggcgctgcgccagcacccttgcgcagc
taccgggtactctgaagaccactccagagcgcacgcccggagcagcaacgggcagaggacattaccaggtataa