



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GOLGA4-golgi autoantigen, golgin subfamily a, 4 Gene

Reference: PTXBC050444

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRPF46; GCP2; GOLG; MU-RMS-40.18; golgin subfamily A member 4; 256 kDa golgin; 72.1 protein; centrosome-related protein F46; golgi autoantigen, golgin subfamily a, 4; golgin-240; golgin-245; protein 72.1; trans-Golgi p230; golgin A4		

For research use only.

GENE DESCRIPTION	golgi autoantigen, golgin subfamily a, 4
GENE ID	2803
ORIGIN SPECIES	Human
ACCESSIONS	BC050444
VECTOR	pDONR223
NCBI SYMBOL	GOLGA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcacgaataactgaagcaaaagatcagcgaggagcagcagcagctccagcaggcgctggctctcgtcaggctcctcaattctcacaacaaagaatgaggagcaggacatcttattacagagcaactgatgaaggtacaccaatagagatcagggtgacacacgtcttgcagacagagctccagctcggctcgccctccgctggagagctttgttcgaagctcgataaagaatctctatcttcctctcttaaaagatcttggatcagaacatctccagagaaatccctgaatcagctgaactcagctgacctgacgtcagacatcttgcagttcctgagttgattgccccccttgataggtatagcgagctcgaagaactctatgaggaattcagactctcacaacaacagatctatcagaaactcttcgaagaagtgaaacgaagcttaagtactcagctcagaggaataattctcagatctccgctcattctatctca

More info: Please contact us for any question: contact@proteogenix.fr