

FLJ30851-hypothetical LOC653140 Gene

Reference: PTXBC050462

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC653140
GENE ID	653140
ORIGIN SPECIES	Human
ACCESSIONS	BC050462
VECTOR	pDONR223
NCBI SYMBOL	FLJ30851
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgccacaaaactgcgagttatgatgaacattcaggccagaaaagttaagagaatggccggagcccgagtcctgttcttaatggaggtattagctagagaagacattgatgaagcagtatgtgcaatattattaaagaaaattccattgtgaaggttacagtcaccattgttgatc
ctctgttcaaaagacagcaagaggtggaagagagagaactggtcttcagtgtagacaggaaaacgacattccataaaagaactgaagaaatagagaagccagctgcatgccagctgccctacttcaacttcaacttcacactgtgtgattccaaaagagtgccataaagcctct
gcaagagccaggagtaaaacttacaatacagtcctgaaaagctcatctatgcagacaagaacagaaaagaaaagagaaaaagacggccgacctaagtcaggctgctgttgaaggcagttccttctcaagctcagccagaagaacaaagtggtgaaagaagggtctgtga
gcagaggcctggggcggtggcatgcagggttgcagcaccacagcagcacatactgttcagaatga