



Reference: PTXBC050526

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein CASC3; BTZ; MLN51; MLN 51; barentsz; cancer susceptibility candidate gene 3 protein; metastatic lymph node 51; metastatic lymph node gene 51 protein; protein barentsz; cancer susceptibility candidate 3		

For research use only.

GENE DESCRIPTION	cancer susceptibility candidate 3
GENE ID	22794
ORIGIN SPECIES	Human
ACCESSIONS	BC050526
VECTOR	pDONR223
NCBI SYMBOL	CASC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]