



Reference: PTXBC050550

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAN6A8; MANA; MANA1; alpha-mannosidase 2C1; alpha mannosidase 6A8B; alpha-D-mannoside mannohydrolase; mannosidase, alpha 6A8; mannosidase, alpha A, cytoplasmic; testicular tissue protein Li 115; mannosidase alpha class 2C member 1		

For research use only.

GENE DESCRIPTION	mannosidase, alpha, class 2C, member 1
GENE ID	4123
ORIGIN SPECIES	Human
ACCESSIONS	BC050550
VECTOR	pDONR223
NCBI SYMBOL	MAN2C1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]