



Reference: PTXBC050553

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BNHS; LNMS; NTE; NTEMND; OMCS; SPG39; iPLA2delta; sws; patatin-like phospholipase domain-containing protein 6; patatin like phospholipase domain containing 6		

For research use only.

GENE DESCRIPTION	patatin-like phospholipase domain containing 6
GENE ID	10908
ORIGIN SPECIES	Human
ACCESSIONS	BC050553
VECTOR	pDONR223
NCBI SYMBOL	PNPLA6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1