

ZNRD1-zinc ribbon domain containing 1 Gene

Reference: PTXBC050608

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HTEX-6; HTEX6; Rpa12; TCTEX6; TEX6; ZR14; hZR14; tctex-6; DNA-directed RNA polymerase I subunit RPA12; RNA polymerase I small specific subunit Rpa12; transcription-associated zinc ribbon protein; zinc ribbon domain containing 1		

For research use only.

GENE DESCRIPTION	zinc ribbon domain containing 1
GENE ID	30834
ORIGIN SPECIES	Human
ACCESSIONS	BC050608
VECTOR	pDONR223
NCBI SYMBOL	ZNRD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgtcatggacctcgccaatacttgcaccagtttcagtcggacctggattctgttcagattgcggctcggtcctctgcccgggctcaggatacggtcacctgtattcgtgtggttcaacatcaacgttcgggacttgaggggaaggtgtgaagactcggttgttcacc
aactggggacagccatgcctatgtcgggtgaggaaggcctgagtgccaggacgtgtggttgacaggcgcctgcctcgtatgtggtcatgaaggaatggcaccacaccagacaaatgcgttcagccgatgaaggcgaactgtcttctacacctgtaccaactgcaagttccaggag
aagggaagactctga