

ASCC3-activating signal cointegrator 1 complex subunit 3 Gene

Reference: PTXBC050681

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ASC1p200; HELIC1; RNAH; activating signal cointegrator 1 complex subunit 3; ASC-1 complex subunit P200; RNA helicase family; helicase, ATP binding 1; trip4 complex subunit p200		

For research use only.

GENE DESCRIPTION	activating signal cointegrator 1 complex subunit 3
GENE ID	10973
ORIGIN SPECIES	Human
ACCESSIONS	BC050681
VECTOR	pDONR223
NCBI SYMBOL	ASCC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctttacctcgtctcacaggagccttgcgttctttcaaatgtcaccaagcaagataattataatgaagaagtggctgacttaagataaagcgaactcatgaacaagttagattgggcctgacatggaagaagataaaaattttgaatgaaaaactggagaagagtaaaatgcaagataaatgaagactaaaagatatattacatgctgcaaacgagatagaagtgaattgtccattccagaagaggaggctggatggaaaaggaggatgagaagatgagcagagcttctgacagattcagaggactaagatga