



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

NUDT10-nudix (nucleoside diphosphate linked moiety X)-type motif 10 Gene

Reference: PTXBC050700

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APS2; DIPP3-alpha; DIPP3a; diphosphoinositol polyphosphate phosphohydrolase 3-alpha; diadenosine 5',5'''-P1,P6-hexaphosphate hydrolase 3-alpha; diadenosine hexaphosphate hydrolase (AMP-forming); nucleoside diphosphate-linked moiety X motif 10; nudix (nucl		

For research use only.

GENE DESCRIPTION	nudix (nucleoside diphosphate linked moiety X)-type motif 10
GENE ID	170685
ORIGIN SPECIES	Human
ACCESSIONS	BC050700
VECTOR	pDONR223
NCBI SYMBOL	NUDT10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr