

FLJ10490-hypothetical protein FLJ10490 Gene

Reference: PTXBC051004

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein FLJ10490
GENE ID	55150
ORIGIN SPECIES	Human
ACCESSIONS	BC051004
VECTOR	pDONR223
NCBI SYMBOL	FLJ10490
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctaaaggttgattcaagcgggggctgtccggaagacgcgtgtgttgaaggtggagtgagcgcccggtggcgaggcgacctcattctgcacccctgcgcgcgcctcgggaactgcacggcgacccccaccgcgactccacgcagacagtagtgcggcc
tgcagggttccccggcgacgaggtaatgttcgtccgccccgccacacagaggcgccactggctccggtgcgtttcagcactctgaggaaaggacttggccttcgcccacagcgtcttaagggtaggcggcggtgtcctcagttctgccccagcactcagacccagac
tgggtccctgcctccgcccctccgctcgactag