

ZNF69-zinc finger protein 69 Gene

Reference: PTXBC051309

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Cos5; hZNF3; zinc finger protein 69; ZNF3		
For research use only.			

GENE DESCRIPTION	zinc finger protein 69
GENE ID	7620
ORIGIN SPECIES	Human
ACCESSIONS	BC051309
VECTOR	pDONR223
NCBI SYMBOL	ZNF69
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctgctgtagtcacaggaggtgtagagaggacccgggacatctgaaagccaggaaatggaggagtggtggcttggctggtatattccagaggaaactctacaaggaagtgtgctggaacttcaggaaactgacctctgtaggaaaaagttggaagaccagaacattgaatat
gagtacccaaacccaggagaacttcaggagtcctagaaaaagaaagtcattgaaattaaagatgacagtcattgtggagaaactttaccagggtccagatgacaggctgaacttcaggagaagaaagcttctctgaaataaatcatgtgacagcttgtgtgtggagaagttggc
ctaggtaactcatcttttaatatgaacatcagaggtgacattggacacaaggcctatgagtatcaggaatatggaccgaagccataa