

FLJ10213-hypothetical protein FLJ10213 Gene

Reference: PTXBC051316

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EBLN-2; endogenous Bornavirus-like nucleoprotein 2; EBLN-1; endogenous Borna-like N element-1; endogenous Borna-like N element-2		

For research use only.

GENE DESCRIPTION	hypothetical protein FLJ10213
GENE ID	55096
ORIGIN SPECIES	Human
ACCESSIONS	BC051316
VECTOR	pDONR223
NCBI SYMBOL	FLJ10213
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtattttcctaattgtatgcttatgttaatttcacagctctttgtttgggtctgtgacagatcttacaagaagatcttttagacatgattcttaacaaataaagaatttaagtcggaaccaattttccacatgtctcatgaagaaggactcacagccagcagccagagatgacgcaa
tggacaggagtggtccctgacctcaaggaagatttgagctatctgggaaaaacagacagatccactggatgcattgggaacccaacccagcattgggataattaaggacattaaaaagcagccaagtctatgtagaccagcacataaatctcatttccacctgtgaccccaagtt
tagtattctgtgttcataattgatgggttacaccaggcattactgagtggtgtgagcaagaggtctaaactgtgttggaatgagaacgaggaagggttactccttatgtagcagattcaagatatgcctaactttattgcccttgagaagtcacagttctccgcaactgtgtgacct
tttgataggcattgcggctggatcaagtataaagttgacacagcagtcaccaagttcagagacgattcaaggcaatgatggcatctattggaagactttcacatggtgagagtgctgctgctaatcagctgcaatgcagaatcagccatagggttgatcagctcaagaccgtgggttga
gaattaatgttcacacttctatttgagactttgaatccctctacacaagctacgaagtcaggttag