



Reference: PTXBC051322

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	volume-regulated anion channel subunit LRRC8A; AGM5; SWELL1; leucine-rich repeat-containing protein 8A; swelling protein 1; leucine rich repeat containing 8 family member A		

For research use only.

GENE DESCRIPTION	leucine rich repeat containing 8 family, member A
GENE ID	56262
ORIGIN SPECIES	Human
ACCESSIONS	BC051322
VECTOR	pDONR223
NCBI SYMBOL	LRRC8A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]