



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

IL9R-interleukin 9 receptor Gene

Reference: PTXBC051337

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD129; IL-9R; interleukin-9 receptor; IL-9 receptor; interleukin 9 receptor		
For research use only.			

GENE DESCRIPTION	interleukin 9 receptor
GENE ID	3581
ORIGIN SPECIES	Human
ACCESSIONS	BC051337
VECTOR	pDONR223
NCBI SYMBOL	IL9R
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatctcgggaagtaactctgctcaagaacgcagacacactctgctcagagaactctgccacggctttcatctgctgtgctgtggtttcaggctgtcacctctcattctagaaaggggccctcagcccagtcctctcagggctgtgacctgtgagagtgaggccctgagcgcagacatggg
cacctgtgctgctgctcagatctcactctgctgtctgtcggagatctctgtcacagggaagcagaagggccaaggtctgaactcctcaactcctcaacaaactcttcagatgctgattgccactgtctgcccagcgtgtgcccagggctgcagccctgctctcttca
caggctctctgcggcacaataatgctacatgtctcgtggggcagtgtagtcacacgtctgtctgccacaggcagcagtgctgtctgcactctgcaatttccatccactctcccaactctcagctctggagagcagcgtctcagctcgtggaccctgcggcagacatcccccagacacagc
aacatcagctcttgcggcactgcatctgcagctggcagctacagctcgtctgtggagccaaatgacacaaactctcagctatgtagctggtccctcaagaagcagggaagggcctgggagcagggccagcagggatcattgtcgggggtgactctggtctatacttgaagcccttgagctggaac
ctgtctctatccatgtagcagctcgtgtctgcatgagctgacccctcagctgaggaattgctgtagagagcaggagcgtctatatacagcagctggagagctggagctgagcagcagcctgtgtctcaggtctcccacagagacaaggccctctgatccacacctgtggtggcagggcaaacacctgtg
ctgtctctctctctctctctcactctcctcagctcccccactctctcaaacctctctcccaactctgatatggggcccccagccggctgtctctcttga