

## OLIG3-oligodendrocyte transcription factor 3 Gene

Reference: PTXBC051352

|                |                                                                                                                                                      |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                  | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                 | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | Bhlhb7; bHLHe20; oligodendrocyte transcription factor 3; class B basic helix-loop-helix protein 7; class E basic helix-loop-helix protein 20; oligo3 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | oligodendrocyte transcription factor 3                                                            |
| <b>GENE ID</b>          | 167826                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC051352                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | OLIG3                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaattctgattcgagctctgtctccagcagagcttcattccggacatggatgagatgtacgtgagggaccaccaccaccaccaccaccaccaccaccaccagagagccgtctcaactcgtctcgtccacgcaggcgcatatgatcagaagatgccggggaaagcctctcggggctgg  
cgccaagggcggggagagagcagcaagtacaaaatcaagaagcagctgtcggagcaggacacacagcttgaggctgaagatcaacggacgcaagcgaagcggatgcacgacctgaacctagccatggacgggctgcgcgaagtcatgccctacgcgcgcatgggctcgg  
tgcgaagctctccaagatcggcaccactcctgtcgcgagaactacatcctcatgtccaccagctccctggaggagatgaagaggctgtgtgagagatctatggggccaccactcggcctttcactgcgggaccgtgggaccactcggcgccaccgccgcacgcggccaactc  
cgtgcacccgggtgcaccccatcttggggggcgctctcatctggcaagcctcgtcaccgctgtccggcgctcactcccgccatcgccaccatcggcctccaccctcgtactcaaggcgccctccacgcgcccgctgcagctgggcagcggtctcagcactgggctggtc  
tgccctgccccctgcacatctgccagatgccggcgccgcacctgtccgtctctccacagccaatggccggctgtcggcgagtcgaaggacttgcgaagtga