

LOC100128510-hypothetical protein LOC100128510 Gene

Reference: PTXBC051366

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC100128510
GENE ID	100128510
ORIGIN SPECIES	Human
ACCESSIONS	BC051366
VECTOR	pDONR223
NCBI SYMBOL	LOC100128510
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgatggtatgctgaaagtggctgctgacttgaccctcagatcagggccttggggagacaaagggcctgcaagagagagctcggtcagagtaaaaaatctctcaaatggactgtctcgatcctaagagctctgctggcacttcgggaattcacctatgatgaagcaggtgga
ccccgtgaggctgtcagcaaaactcaagaattatgcatctatggctgaagccagagatccactcaaaagagcagatctggaactgctgtgctggagcagttcctgactattctgccaggagacagaccagatgcagaagccatccacagacattgaggaggctgtgct
ctggtagaacacttcagaggggaatctggtcaaacatggaatgggtgagaagaagattcctgacatgtacagtgaatatggatgcaggtggatag