

DKFZp434K191-hypothetical locus DKFZp434K191 Gene

Reference: PTXBC051721

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical locus DKFZp434K191
GENE ID	29797
ORIGIN SPECIES	Human
ACCESSIONS	BC051721
VECTOR	pDONR223
NCBI SYMBOL	DKFZp434K191
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaaagacacagagaaggatgtgacatgcctggccatggagcactctgagatctcatcgtggacaccactgccacacccatcccgtctgcgcaggccgacactcactgacgttgaaggctgccttcagtcctggatgtctgcggccacccacagatgcggtagatgc
ccacctctccatgcctcggcgtcgtatctctccacgcactggcgacgatgaggcacttggacctctctcctcggggaggagggaatgttctcagtgctcctaacagccctgttggccataa