

## SLC44A5-solute carrier family 44, member 5 Gene

Reference: PTXBC051740

|                |                                                                             |                     |                |
|----------------|-----------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CTL5; choline transporter-like protein 5; solute carrier family 44 member 5 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | solute carrier family 44, member 5                                                               |
| <b>GENE ID</b>          | 204962                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC051740                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | SLC44A5                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgatcttctactgtgtatttggctacattgttttagaacttggcctgggtacatggggacccagaagagcagcctctacagacagccaggccacttttggccagaaggcactcccaatgagaacaagaccattttgtttactttaacctgttacgctgtaccagtcctccgt  
attgctaaccctacagtgcctaccacacagatctgtgtcctcaagtggccagaaaaattttaacctatgtggaatgcaactttgtacacaaaagacaaaagctactgggaagactaccgctcagttctgtaagaccactgtaagcctgtgaagctctcacacagctttactggatgatgatt  
gtccaacagcgattttccacgaaaccttttccagagatgttccctgactctctacaaaaatggcactttaacaatagggaatgaagatgatgtttcaagatggaaatggaggacaaagtggtgtaaacctcgggattgctgcaaatggtatcaataaactcttgatgcaaaagtcactt  
ggattgaaagtgttgaaactatgcaagaacttggatgtgattctcattggcctgacgattgccatggcttagttggatatttttgatacttctgaggttcatagtcggatgcctcttctgggtcttcatgattgggtgtgattggaattataggttatggaatatggcactgttaccagcagtacacc  
aatcttcagggaacgccaaagtctgtatttaactatctatgacatcgggattcagactaataaagcatgtacttgaactgcaacaacatgtgtcacattatgataaactctgcatcattgaagtgaattgcatctcctcatgctgatcttctcagggaatcgaatccgagtcgccattatcctgctgaa  
ggaagggaagcaagccattgatatgttctagtagtctatccagctttaacttttcttctcaatctgatttgcactgggtcgtgacagcagttttcttggcgacatcgggggtacctgtatacaaaatcatagctccagggggcattgtatacatgaaaatcaaacctgtgaccc  
agagatttttaatacaactgaaatggccaaagcttgcctggggctctgtgtaacttgccttctatggtggaagagcctgttaccatcagtagatccctaccctccatgatacaactatttcttctctggcctataaactctgctcattgtagtcagtcgcgcctgtctggtgactgctactt  
attactgggccaatgaaaaacctgatgacatccacacatccacttttactgacttggagcagcctacagatcacacagcagcctcctgacttggatcttaattatgcaattcaaatgtttaaattgtactagaatacttgaccaccgtcttaacgtaccagaacacattgtctaaa  
ttcctacaatgctgctgagatgctgcttctgtgtttggaaaatgcaataaagttttaaacagaaatgcctatattatgattgcaatataatggcagaaactctgcaggcagcaaaagatgcttcaatctgctgatgagaatgtttgaaagtgcagttacagatgaagtacatacttattat  
tctggggaaactctagttgctggaagtataggtgttctggcctctctattctcacaaaagactgccaagtgtgacacagaccagcactttaaattactactgggtaccttctgacagtcattttgggtcttacctgattgcacatgggttcttcagcgtctatgcaatgtgttgaaca  
atttcatctgcttctggaagatttagaaagaatgatgttctactcgaagacctattatgtgagtcacaccttctggaagatttccagaggagaaatccacaactaggaagcagtag