



Reference: PTXBC051743

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZFP91 zinc finger protein; zinc finger protein homologous to Zfp91 in mouse; E3 ubiquitin-protein ligase ZFP91; DMS-8; DSM-8; FKSG11; PZF; ZFP-91; ZNF757; drug resistance-associated sequence in melanoma 8; penta zinc finger protein; zinc finger protein 75		

GENE DESCRIPTION	zinc finger protein 91 homolog (mouse)
GENE ID	80829
ORIGIN SPECIES	Human
ACCESSIONS	BC051743
VECTOR	pDONR223
NCBI SYMBOL	ZFP91
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]