



Reference: PTXBC051744

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACCPN; KCC3A; KCC3B; solute carrier family 12 member 6; K-Cl cotransporter 3; electroneutral potassium-chloride cotransporter 3; potassium chloride cotransporter 3; potassium chloride cotransporter KCC3a-S3; potassium-chloride transporter-3a; potassium-ch		

For research use only.

GENE DESCRIPTION	solute carrier family 12 (potassium/chloride transporters), member 6
GENE ID	9990
ORIGIN SPECIES	Human
ACCESSIONS	BC051744
VECTOR	pDONR223
NCBI SYMBOL	SLC12A6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]