

WDR54-WD repeat domain 54 Gene

Reference: PTXBC051753

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WD repeat-containing protein 54; WD repeat domain 54		

For research use only.

GENE DESCRIPTION	WD repeat domain 54
GENE ID	84058
ORIGIN SPECIES	Human
ACCESSIONS	BC051753
VECTOR	pDONR223
NCBI SYMBOL	WDR54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcgctgggagcgtccattccccctgcgaggctcggcgcgcctgtgcaacaacctcagtgctgcagctgccggctcgaacctcagttttggcgtggttcattggaccaagcggccagcttcagcgtgctcctgagggtgtgcccttgccccagcggcagctccacg
ctaaggagggtgctggagtgagtcctcccaattatcactcaggtgaggcatggagctgggagtgatgttgctgaacccagtcctacccccagcctttctctccaaaccttcaggagcaggcatgtcctaacctcagaatctcccagcatgtcctaa