



Reference: PTXBC051765

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CVID12; EBP-1; KBF1; NF-kB1; NF-kappa-B; NF-kappaB; NFKB-p105; NFKB-p50; NFkappaB; p105; nuclear factor NF-kappa-B p105 subunit; DNA-binding factor KBF1; NF-kappabeta; nuclear factor NF-kappa-B p50 subunit; nuclear factor kappa-B DNA binding subunit; nucl		

For research use only.

GENE DESCRIPTION	nuclear factor of kappa light polypeptide gene enhancer in B-cells 1
GENE ID	4790
ORIGIN SPECIES	Human
ACCESSIONS	BC051765
VECTOR	pDONR223
NCBI SYMBOL	NFKB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]