

STX17-syntaxin 17 Gene

Reference: PTXBC051790

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	syntaxin-17; syntaxin 17		
For research use only.			

GENE DESCRIPTION	syntaxin 17
GENE ID	55014
ORIGIN SPECIES	Human
ACCESSIONS	BC051790
VECTOR	pDONR223
NCBI SYMBOL	STX17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgaagatgaagaaaaagtgaattacgccgtcttgaaccagctatccagaattcattaaagatagtaatccaacagacctggaaggttaagaaagcaccagataaatattgagaagtatcaaaggtgcagaatctgggacaagttgcatgaagagcatatcaatgcaggacgtac
agttcagcaactccgatccaatccgagaaattgagaaactttgttgaagtcggaaggatgacctagctactctgaagagaatgatagatcctgttaagaagaagcatcagcagcaacagcagaattctccaactccatttgaatctgtagaagaacttaagaagcaatttaagatga
agaaactttgctacagcctcctttgaccagatccatgactgttggtggtatgtattgtgctaa