



Reference: PTXBC051816

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM113A; C20orf81; bA12M19.1; PC-esterase domain-containing protein 1A; family with sequence similarity 113, member A; sarcoma antigen NY-SAR-23; PC-esterase domain containing 1A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 113, member A
GENE ID	64773
ORIGIN SPECIES	Human
ACCESSIONS	BC051816
VECTOR	pDONR223
NCBI SYMBOL	FAM113A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]