

PCGF5-polycomb group ring finger 5 Gene

Reference: PTXBC051845

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNF159; polycomb group RING finger protein 5; RING finger protein 159; ring finger protein (C3HC4 type) 159; polycomb group ring finger 5		

For research use only.

GENE DESCRIPTION	polycomb group ring finger 5
GENE ID	84333
ORIGIN SPECIES	Human
ACCESSIONS	BC051845
VECTOR	pDONR223
NCBI SYMBOL	PCGF5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctacccaaggaacacttgggaaagattttaacccatcattacgtctatatctgtaaagggtatctgatcaagccaacacagtgacggaatgcctccatactctgtaagacttgattgttcagcactttgaagatagcaatgattgcccaagggtggcaaccaagttcatgagac
aaatccattagaatgttgaggttgacaatacattagaggaaattatattaaagctggcctcgactacgagaacaagaacttgagcgtgaatctgaattttggaagaaaaataagcctcaagaaatggacaagatgatactcaaaagctgacaaaccgaaatgagatgaagaagggtgat
gaaaatgaagatgataaagattatcacagaagtgacccacaattgctatctgtctagattgtttacgaataatgggcaatcaggggacaatgtagtaaaagggtttaatgaagaaattcattcgatgttctacacgtgtgaactgtgggaactattaaaaaatttcaagtttaaaactaaaactcca
agttcttatgagttggatgtctgtgcaatggtgaaattatggggaaggatcactatggaattcatctacatgacaagatggcgactaagaggcgaactttcggtgtctgaactgctcagcttcgaagctgctctcaggatggccctttgtatcagtcataccctatggtactgcagtatcg
accaagaattgattcgggttag