



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

WDR54-WD repeat domain 54 Gene

Reference: PTXBC051848

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

WD repeat-containing protein 54; WD repeat domain 54

For research use only.

GENE DESCRIPTION	WD repeat domain 54
GENE ID	84058
ORIGIN SPECIES	Human
ACCESSIONS	BC051848
VECTOR	pDONR223
NCBI SYMBOL	WDR54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr