



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RPAIN-RPA interacting protein Gene

Reference: PTXBC051849

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HRIP; RIP; RPA-interacting protein; RAP interaction protein; nuclear transporter; RPA interacting protein

For research use only.

GENE DESCRIPTION	RPA interacting protein
GENE ID	84268
ORIGIN SPECIES	Human
ACCESSIONS	BC051849
VECTOR	pDONR223
NCBI SYMBOL	RPAIN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atfgcggagctcgttgaggtctccgccgccctccctgtacaactgggtggctgccgcctfggaagaggcttccggcagagatgcctggagagaatgagaacacgccgggacaggctctaacaaggctaccgcaggctggaacagctgggccagggaattctcagaacagcttt
ctagtccaagaggtgatggaagaagatgggaatgcttgcagctcagctggagaatgcttcagaagacttggctcagttggggagagctgatagacatggctgtgctggaggaaattcaacaggagctgatacagaagaagcagctccatcatcagcagatgatgagaagagcttgcagtgtgatg
aaaagctgtcacatcatctgctgcagtgaggctgaacacacactcatctctgtatgacaaatgacaaacacggctggtggtgtgctgatacagctgcttcacatcccatctattctctgtagacagcagaagcagctgctgctgtttagagggtagtataaatg
acacagctgtcactatccccacacactgaatttccatcactcgcgaagaacacagatggacagctctctcatatgcgtctctgcctctgtacatctggctctgctcctgag