



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C19orf46-chromosome 19 open reading frame 46 Gene

Reference: PTXBC052573

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf46; DFNB76; Nesp4; nesprin-4; deafness, autosomal recessive 76; nuclear envelope spectrin repeat protein 4; spectrin repeat containing nuclear envelope family member 4		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 46
GENE ID	163183
ORIGIN SPECIES	Human
ACCESSIONS	BC052573
VECTOR	pDONR223
NCBI SYMBOL	C19orf46
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr