

NCR3-natural cytotoxicity triggering receptor 3 Gene

Reference: PTXBC052582

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1C7; CD337; LY117; MALS; NKp30; natural cytotoxicity triggering receptor 3; NK-p30; activating NK-A1 receptor; activating natural killer receptor p30; lymphocyte antigen 117; natural killer cell p30-related protein		

For research use only.

GENE DESCRIPTION	natural cytotoxicity triggering receptor 3
GENE ID	259197
ORIGIN SPECIES	Human
ACCESSIONS	BC052582
VECTOR	pDONR223
NCBI SYMBOL	NCR3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctggatgctgttgcctccttgatcatgggtccatccagatcctgtgctctctgggtgtccagccccctgagattcgtacctggaaggatcctctgcttctgcccctgtccttcaatgccagccaaggagactggccattggctccgtcacgtggtccgagatgaggtgttcca
gggaaggaggtgaggaatggaacccagagttcagggccgctggccccactgtcttctccctgttctccatgaccaccaggtgagctgcacatccgggacgtgagggccatgacgccagcatctacgtgtgagagtgaggtgtggtgctggccttggtgtcgggacagggaatg
ggactcggctgtgtgggagaaagaacatcctcagctagggtgtgtacagtcctcctcctcgggtggattctatgctgtcagcttctctctgtggccgtgggcagcacgtctattaccagggcaaatgtctgacctggaaagggtccaagaaggcagctgccggtgtgtccagcg
cccccccaccacatgtgggagctcagcacatctgttccccagtcaccaggagctga