



Reference: PTXBC052590

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A2HS; AHS; FETUA; HSGA; alpha-2-HS-glycoprotein; alpha-2-Z-globulin; ba-alpha-2-glycoprotein; fetuin-A; alpha 2-HS glycoprotein		

For research use only.

GENE DESCRIPTION	alpha-2-HS-glycoprotein
GENE ID	197
ORIGIN SPECIES	Human
ACCESSIONS	BC052590
VECTOR	pDONR223
NCBI SYMBOL	AHSG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1