

CMC1-COX assembly mitochondrial protein homolog (S. cerevisiae) Gene

Reference: PTXBC052644

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf68; COX assembly mitochondrial protein homolog; C-x(9)-C motif containing 1; COX assembly mitochondrial protein 1 homolog; CX9C mitochondrial protein required for full expression of COX 1; cmc1p; mitochondrial metallochaperone-like protein; C-X9-C mo		

For research use only.

GENE DESCRIPTION	COX assembly mitochondrial protein homolog (S. cerevisiae)
GENE ID	152100
ORIGIN SPECIES	Human
ACCESSIONS	BC052644
VECTOR	pDONR223
NCBI SYMBOL	CMC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctcgacccgcagaccagcatctcagacatgtcgaataatgtttgatccctaaataatgagagaaaaggccaagagagggtgtctgaacaagtccaagatttaccatgttgcaagaactctggagttcttatggtagtaaaatgccggaagaaaattctgcattgaaag
aatgtctaactgcttactataatgatccagcctttatgaagaatgcaaatggaatacctgaaggaaaagggaagaattcagaaaaactggaattcctacaagaaaaggctacagaagcttccaacagcatgtag