

MRPL55-mitochondrial ribosomal protein L55 Gene

Reference: PTXBC052806

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AAVG5835; L55nt; MRP-L55; PRO19675; 39S ribosomal protein L55, mitochondrial; L55mt; mitochondrial ribosomal protein L55		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L55
GENE ID	128308
ORIGIN SPECIES	Human
ACCESSIONS	BC052806
VECTOR	pDONR223
NCBI SYMBOL	MRPL55
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgtgggcagctgttggcggctgaggcagagcaccgtgaaggccacggacgtgcactccgcgcctgcacacatcctctggcgagctgacagcagggcctcactcactcgtgtgcaccgccaggttatgcagactctacccgtgctgctggtgaagcag
gatggctccaccatccacatccgctacagggagccacggcgatgctggcgatgccatagatctggacaccctgtctctgaggagcgccggcaggtgctggaagcgtgaggctcagctccagtcgaggaaggagtacgagcaggagctcagtgatgactgcatgtggagcg
ctaccgacagttctggaccagaccaagaagtga